



SOP
TA Instruments SDT-Q600

I-02	
Revision	03
By	Sami Vuori
Issued	22/4/2022

Standard Operating Procedure
TA Instruments SDT-Q600 TGA/DSC

Table of Contents

1. Introduction.....	2
3. Measurement.....	2
4. Data export.....	9
5. Performance data.....	10
6. Version log.....	10



SOP

TA Instruments SDT-Q600

I-02

Revision	03
By	Sami Vuori
Issued	22/4/2022

1. Introduction

The thermal analysis studies the relationship between a sample property and its temperature as the sample is heated or cooled in a controlled manner. The thermogravimetric analysis (TGA) evaluates changes in the sample mass as a function of temperature, while the differential scanning calorimetry (DSC) measures the differences in the heat flows between a sample and the furnace as well as a non-reacting reference material, and the furnace as a function of temperature. The TGA-DSC measurements are often used to obtain information about the physical and chemical phenomena which may be associated to mass changes and /or heat changes. This equipment provides simultaneous TGA-DSC measurements from ambient temperature to 1500 °C.

2. Safety

Normal laboratory precautions: a lab coat, gloves, and goggles.

3. Measurement

3.1. Switch on the computer and the equipment from the switch at the back of the instrument.



3.2. Log in to the computer (username **kemlaite**, password **Maja123paja**) and right click the program **TA Instrument Explorer shortcut** → *Run as administrator*.



SOP TA Instruments SDT-Q600

I-02

Revision


03

By

Sami Vuori

Issued


22/4/2022



3.3. When the machine has initialized (takes a few minutes), double click the instrument in the window; there is a red circle on the icon until it's ready.



You may see some errors; just click **OK** on them.



3.4. Decide which gas you will use in the analysis and open the gas flow accordingly: use only the valve marked in the photo with red arrows; do not touch the other two marked with the red **DO NOT TOUCH** labels!



SOP

TA Instruments SDT-Q600

I-02

Revision	03
By	Sami Vuori
Issued	22/4/2022



Gas #1 = N₂ (Typpi N₂), **Gas #2** = pressurized air (paineilma).

3.5. In the program window, click **Control** → **Gas** → **Gas 2** (or **Gas 1** if you're using nitrogen).

3.6. In the **Summary** tab, change the data:

Summary | Procedure | Notes

Procedure

Mode: SDT Standard

Test: Custom

Sample Information

Sample Name: Sample name

Pan Type: Alumina

Comments: Air

Data File Name: \\Ty1506087\ta\Data\SDT\Sami\hackmanite.001

Network Drive

Mode: SDT standard

Test: Custom (you may choose some other if you want, but this is generally used)



SOP

TA Instruments SDT-Q600

I-02

Revision

03

By

Sami Vuori

Issued

22/4/2022


Sample name: [Your sample name]

Pan type: Alumina


Comments: Air (or “N2 gas” if using nitrogen)

Data file name: make a file for this measurement by clicking the icon on the right side

3.7. Go to the **Procedure** tab and check the parameters:



If you want to change the heating parameters, click the **Editor** button and change the values by double clicking them.



3.8. Go to the **Notes** screen. Choose the gas to be used in the measurements and check that the flow rate is set to **100 ml/min**.



SOP

TA Instruments SDT-Q600

I-02

Revision

03

By

Sami Vuori

Issued

22/4/2022

The screenshot shows the software interface with the following sections:

- Notes:** Operator: Sami; Extended Text: (empty)
- Mass Flow Control Settings:** Sample: #2 - Air; Flow Rate: 100 mL/min
- Auto Analyze:** Autoanalyze; Analysis Macro: (empty)

Check that the gas flow (**Sample purge flow**) is actually 100 ml/min in the upper right side of the window.

Signal	Value
Temperature	18.05 °C
Weight	-183.538 mg
Weight percent.	100.00 %
Reference Weight	-187.513 mg
Heat Flow	-3.620 mW
Temp Difference (°C)	-0.42 °C
Temp Difference (µV)	-2.392 µV
Heater Power	0.00 W
Set Point Temp	0.00 °C
Sample Purge Flow	99.97 mL/min

3.9. Open the furnace by clicking **Control** → **Furnace** → **Open**.



3.10. **CAREFULLY** place two crucibles onto the arms (they break easily and are very expensive!).




SOP


TA Instruments SDT-Q600

I-02


Revision	03
By	Sami Vuori
Issued	22/4/2022



3.11. Close the furnace by clicking **Control** → **Furnace** → **Close** and wait for it to close.



3.12. Tare the balance by clicking the scales icon.



3.13. Now open the furnace (3.9.) again and add ~10 mg of your sample to the sample crucible (leave the reference cup empty). In the right part of the window you can check how much sample you have added. Close the furnace (3.11.) after you have enough. After this, in case you need the heat flow curve, wait until the heat flow has stabilized (after ~1 h). If you only need the mass loss data, you can start the analysis when the mass has stabilized.

3.14. Start the measurement by clicking the green start button.




SOP

TA Instruments SDT-Q600

I-02


Revision	03
By	Sami Vuori
Issued	22/4/2022




3.15. Write your run details to the instrument's notebook.

3.16. When the analysis has finished (wait until the furnace temperature is below 50 °C), open the furnace (3.9.), remove the crucibles and close the furnace (3.11.).

3.17. Shut down the instrument by clicking **Control** → **Shutdown Instrument**.



Click **Start** in the next window.



3.18. When the instrument's display shows "*Shutdown complete, it is now safe to reboot your computer*", turn it off from the back (3.1.)

3.19. Turn off the gas flow from the valve you opened in 3.4.

3.20. Clean your crucibles by scraping loose material off, submerging them in aqua regia and finally firing the crucibles at 1500 °C.



SOP

TA Instruments SDT-Q600

I-02

Revision

03

By


Sami Vuori

Issued


22/4/2022

4. Data export

4.1. Double click the TA Universal Analysis program on the desktop.



4.2. Select your file by clicking **File** → **Open** and double click your file in the **C:\TA\SDT\Data** folder. Click **OK** in the next window.





SOP

TA Instruments SDT-Q600

I-02

Revision

03


By

Sami Vuori


Issued

22/4/2022

4.3. Export the data by clicking **File** → **Export Data File** → **File and Plot Signals**.



4.4. Tick/untick the following items and click **Finish**, then save it as a text file to your folder.



5. Performance data

Balance sensitivity: 0.1 µg

TGA sensitivity: 0.001 °C

Calorimetric accuracy/precision: ±2 %

6. Version log

Revision	By	Date	Notes
01	Hellen Silva Santos	3/14/2014	Initial version
02	Hellen Silva Santos	22/03/2018	
03	Sami Vuori	22/04/2022	Updated the SOP for Aurum.